

ReadyWrap® - The easy solution for self-care compression



Easy to use³, with colour-coded VELCRO® fasteners



- Front-fastening short straps aid patient application and readjustment^{3,4}
- Overlapping straps minimise user application error and can prevent undue loss of pressure if the patient does not apply the product perfectly¹
- Low-profile design, allowing patients to wear their usual footwear²

Step 1. Apply ReadyWrap® Liner



Step 2. Apply foot piece



Step 3. Apply calf piece



ReadyWrap® Extender Strap



All products are available in black or beige

Helpful Hints

For comfort and protection you can wear a liner underneath your ReadyWrap®. ReadyWrap® is supplied together with a suitable liner for your convenience, however, any non compressive sock or liner may be used.

You may prefer to apply your lower extremity ReadyWrap® garments while 'long' sitting (with legs stretched out in front). However, some find it easier to put on ReadyWrap® whilst sitting in a chair or standing. We encourage you to see **what method works best for you**.

The full ReadyWrap® range includes, Thigh, Knee, Calf, Foot and Toe pieces as well as the Liner.

All ReadyWrap® garments can be worn in the day and/or night. Please follow your clinicians guidance.



A full holistic and vascular assessment should be carried out before treatment. Please refer to local and national guidelines.

The ReadyWrap® - How to Apply video is also available to view on the L&R Medical UK YouTube channel. Please scan the QR code below ▼

For further information, please contact your local **L&R representative** or **L&R Customer Services Team** on **08450 606 707**

