

Cellona® Shoe Sizer

Place the bandaged foot on this guide and select the closest size.
NOTE: If the foot is between 2 sizes then choose the larger size.

To prevent tripping, use scissors to trim the hook-and-loop ankle fastening tape, as required.

TOE
.....▲ XL

L

M

S

XS

P

HEEL
.....▲

M1696 V2.1



People.Health.Care.

For more information call **08450 606707**
to speak directly with a Customer Care Advisor
or visit **www.Lohmann-Rauscher.co.uk**

| Cellona® Shoe Size | Bandaged Foot Circumference | PIP Code |
|-----------------------|--------------------------------|-------------|
| XL | up to 40cm | 385-9501 |
| L | up to 40cm | 385-9493 |
| M | up to 40cm | 385-9485 |
| S | up to 34cm | 385-9477 |
| XS | up to 28cm | 385-9451 |
| P | up to 25cm | 385-9469 |